

INSTRUCTIONS FOR USE


The Translock™ is a durable, lightweight, and easy-to-use latching device offering cross functionality for wheelchair user transfers. The Translock™ can either be attached to a manual wheelchair or to the POST as shown in the instructions.

The Translock™ can be operated from the wheelchair by simply docking the Translock to one of the companion attachments. With the locking device closed this will prevent the wheelchair from tipping or moving out from under the chair user during transfer.

After returning to the wheelchair release the Translock simply by pushing on the wing tabs to open the locking mechanism. Once latched open the Translock™ may be closed by pressing the release button opposite of the wing tabs. The Translock can then be returned to a stowed position by rotating the Translock back against the wheelchair frame.

If used on the POST configuration just pull the wheelchair up to the Translock™ and make sure its securely locks and attaches to the wheelchair.

Always make sure you lock the wheelchair brakes before transferring.



LIMITED WARRANTY

This warranty covers any defects or malfunction of your new Translock™ device and family of products. The limited warranty is good for one year from the date of your purchase and terminates if you sell or otherwise transfer the product(s).

Transfer Solutions, LLC will replace any device that fails and is returned to the manufacturer for inspection. If the product is deemed to be a failure and not caused by misuse, customers will receive a replacement device at no charge plus shipping and handling.

To be eligible for this warranty return your device with a copy of your proof of purchase and written description of the problem to:

Transfer Solutions, LLC
 311 East Main St.
 Suite 417
 Galesburg, IL 61401

This warranty does not cover product uses of other than intended in product instructions.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

INSTRUCTIONS
 A family of products
 manufactured by:

www.thetranslock.com
 Patent Pending