

Outdoor Stairlift



WARRANTY

- 1YR - Battery
- 2YR - Parts
- 3YR - Drive Train & Circuit Board

- Lowest priced outdoor unit
- Plated rack & pinion gear
- Anodized aluminum rail
- Protective coating on circuit board
- Plated track brackets
- Plated chassis frame
- Stainless steel hardware
- Heavy duty canvas cover
- Sealed switches & waterproof connections
- Outdoor grade powder coating
- Neutral gray color



- Folds to within 11" of the wall
- Attaches to stairs, not the wall
- Narrow, low profile track (2.5" tall x 6.5" wide)
- Obstruction sensors on chassis & footrest with audio beep notification
- Battery power allows 20 trips during power outages
- Top & bottom landing charge stations
- Audio beep notification when unit is not seated on charge stations
- Rocker switch conveniently located on either side of seat
- Smooth starts & soft stops

- Constant pressure controls
- Wireless remote controls for upper & lower landings
- Chassis mounted install switch simplifies installation
- Quickly installs on left or right side of stairway
- Rack and pinion drive with concealed gear rack
- Low-to-the-floor footrest
- 3" to 6" horizontal extension at top landing
- 25 to 45 degree angle applications
- 350 lbs capacity

