




## MobileCharge24EX

---

### MobileCharge24EX (XLR Connector) Connector Compatibility

Compatible with most Power Chairs & Mobility Scooters:

Pride, Hoveround, Permobil, Drive, Active Care, CTM, Sunrise, Heartway, Wu's Tech, Bruno, Merits (all 2006 and later models), Golden (except "Companion"), Shoprider (except Sunrunner Deluxe and Lithium-ion battery models), Invacare (except Lynx L3), Ranger (except Spirit), Quickie (OK for models with off-board charger), PaceSaver (works with Atlas Plus 3, Eclipse, Esprit). May not be compatible with some EV Rider, Zip'r, E&J, or Worldwide Seating products. Works with most other Mobility Scooters and Power Chairs with off-board chargers & XLR connector.



Mobility Electronics, Inc.  
2419 Irma Way  
Castro Valley, CA 94546  
(510) 317-1775

[sales@chargeone.com](mailto:sales@chargeone.com)



Assembled in the USA





Check that your Power Chair or Mobility Scooter has the circular 3-pin “XLR” type connector.

## YOUR 3-PIN XLR CONNECTOR LOOKS SIMILAR TO THIS:



Cable Plug  
(End View)

Battery Socket  
(Opening)



Pin 1 (+)  
Pin 2, 3 (-)