



# KM-5000TP

*Reclining, high strength, lightweight aluminum detachable-arm wheelchairs*

- Folding 6061 T-6 Aircraft-grade aluminum frame in black powder coat finish
- Medicare Code: K0004+E1226+E0966
- Backrest reclines from 90°~163° with anti-tippers
- Padded detachable full length armrests with side panels
- 12” rear polyurethane tires
- Flat free front casters
- Adjustable toggle wheel locks.
- High-grade flame retardant, breathable nylon upholstery
- Limited lifetime warranty on frames and crossbraces, one year on parts\* .
- Frame Color: Black

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open (in.)	weight riggings (lb.)	weight capacit (lb.)
KM-5000TP	18	18	19.5	33 (18.5+14.5)	50.5	27	33	250

Model No.	Description
KM-5000TP	18”, aluminum, 180° reclining back, detachable full length armrests, elevating legrests, anti-tippers, adj. height pillow



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Tel.626-581-2235 Fax.626-581-2335