



KN-880 Series

*reclining, lightweight,
detachable arm wheelchairs*

- * Medicare Code: K0003+K0195+E1226+E096
- * Folding chrome plated steel frame
- * 180° full reclining back with anti-tippers
- * Padded detachable full armrests
- * Detachable elevating legrests
- * 24" x 1" mag wheels
- * 8" x 1" caster wheels with rubber tires
- * Adjustable length leg supports & footplate
- * High-grade flame retardant nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KN-880	18	16	20	29 (+)	49	25 ½	50	250
KN-880N	16	16	20	29 (+)	49	25 ½	50	250

Model No.	Description
KN-880	18", 180° reclining back, detachable full arms, elevating legrests
KN-880N	16", 180° reclining back, detachable full arms, elevating legrests

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Seat Belt • Leg Straps • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335