



## KM-802F

*aluminum high strength, lightweight  
flip-back-arm wheelchairs*

- Medicare code k0004
- 6061 T-6 Aircraft-grade aluminum frame
- Black powder coat finish
- Padded flip-back armrests with side panels
- Swing-away detachable footrests with heel loop
- 24" x 1" wheels with high profile polyurethane tires
- 8" x 1" flat free caster wheels
- Adjustable length anti-slip footplates.
- High-grade flame retardant, breathable nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight Capacity (lb.)
KM-802F18	18	17	19	17	37	25 3/4	29	255
KM-802F16	16	17	19	17	37	23 3/4	29	255

Model No.	Description
KM-802F18	18", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
KM-802F16	16", aluminum, high strength lightweight, flip-back-arms, swing-away footrests



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335