



MVP-502

Ergonomic V-Seat Recliner

- Patented Ergonomic/Anti-Sliding “V-Seating” prevents sliding, sores and relieves overall pressure on your body!
- Backrest reclines smoothly from 93-160 degrees
- Rigidifying/Folding push-bar provides and strengthens back support →



- Ergonomic Armrest allows the user to adjust the armrest position when reclined →



- “Traveling-Axle” Leg-rest automatically stretches smoothly to reduce pressure/compression on abdomen.
- 3D Elevating Leg-rest support Ergonomically Designed to allow greater micro-adjustment: up/down, left/right or inward/outward →



- Extending handbrakes equipped!
- Strap-Tension adjustable backrest & seat for ultimate comfort adjustments and stabilization for each user!
- Anti-Microbial/Ergonomic Upholstery (Seat, Backrest, Headrest & Paint) Made from Aegis Microbe Shield; Also pressure relieving 3cm thick breathable material, offers superb comfort throughout the body
- Folds Ultra Compact and Weighs only ~ 44lbs w/ all fixtures →



items	F/R wheels	Seat	Backrest	Headrest	Armrest	Handbrake	Leg rest	Footrest	Push bar	Backrest	Seat	Max. User Weight
MVP-502	16"	7"/14"	16"	16"	23"	14"	9-11 1/2"	45 1/2"	20"	21"	37 1/2"	33 lbs 250lbs
MVP-502	18"	7"/14"	18"	18"	25 1/2"	14"	9-11 1/2"	48"	20"	21"	37 1/2"	33 lbs 250lbs



KARMAN HEALTHCARE Inc.
 19255 San Jose Avenue. City of Industry. CA 91748-1418
 Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335