


# Koala Klip<sup>®</sup>

**T**he Koala Klip, is a flexible, circular clip that organizes up to six medical tubes, cables and/or IV lines. It consists of a primary clip that can be secured onto IV poles, bedrails, ventilator tubing or used free-floating.

Koala Klip enables medical personnel to keep tubes, cables and IV lines tangle-free. In one quick visual glance, staff can recognize one tube from another. The outer clips are varied in size to accommodate numerous diameters of IV lines, respiratory lines and circuit lines. The Koala Klip is a valuable time saving clip for medical personnel, providing organization and improved patient safety.



Patents Issued D 508, 842

## Koala Klip Features and Benefits

- Organizes up to (6) IV, respiratory and/or circuit lines
- Can be secured to IV poles, bed rails, ventilator tubes or used free floating
- Decreases the risk of mismedications
- Lifts lines, cables and tubes off the floor diminishing the risk of contamination
- Reduces tubing and line entanglement and the need to rework lines
- Minimizes patient entanglement and line impingement
- Increases nursing operational efficiency minimizing medical errors
- Reduces the risk of injury to patients and staff due to tripping
- Non-Sterile and Latex Free
- Disposable, versatile & economical

The Koala Klip is supplied bulk packaged and non-sterile

Item #	Description	Qty/Box
KK1000	Koala Klip Disposable IV Clip	100



Advent Medical Technologies  
P.O. Box 702433  
Tulsa, OK 74170-2433  
Tel: (866) 690-2665 Fax: (918) 362-5003  
[www.adventmedicaltechnologies.com](http://www.adventmedicaltechnologies.com)