

# IMPROVED! QUIETER!

# The most sophisticated heat and massage system is now available from Golden!

## Experience this luxurious new massage system featuring:

- Dozens of massage settings covering three different massage zones: shoulders, lumbar and seat
- Four massage motors provide relief to stiff, tired and strained muscles
- 15 and 30 minute massage timers
- Activate any or all of the three massage zones while in Wave or Pulse massage settings
- Independent adjustable massage intensity and speed settings

## Relax with soothing low and high penetrating heat

- Independent 30 minute heat timer
- Soothe back muscles with our warming heat system
- Low and high heat settings

Press PWR to turn On or Off.

15 Minute Timer:  
Single Press (GREEN)

30 Minute Timer:  
Press & Hold (ORANGE)

Press LO or HI to turn  
HEAT On or Off.  
Independent 30  
Minute Timer

Press any or all ZONE keys  
to turn any Zone On or Off.

Press SELECT, WAVE or PULSE  
for desired operating mode.

Entire display is backlit for  
easy use in low light

Press - (Left) to decrease  
or + (Right) to increase  
INTENSITY and SPEED.



Lead times and upcharges will apply.

Heat and massage system not available on Capri PR200.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518  
Tel: (800) 624-6374 • Fax: (800) 628-5165 • [www.goldentech.com](http://www.goldentech.com)

SHT-DHV.081213