

AmeriGlide 111TS Travel Scooter



The AmeriGlide Travel Scooter's lightweight, compact design makes it convenient and easy to take along with you. Standard features like scratch-resistant plastic, an easily adjustable tiller and off-board battery charger make the 111TS fun and simple to use.

**Travel
Scooter**

Specifications

Model:	AG111TS
Weight Capacity:	250 lbs.
Total Weight:	84 lbs.
Maximum Speed:	Up to 4.25 mph
Ground Clearance:	1.5"
Turning Radius:	32"
Overall Length & Width:	37.75" x 18.625"
Seat to Ground:	19.75" - 17.75"
Seat Dimensions:	Width: 16.5" Depth: 15" Height: 12"
Rear Anti-tips:	Standard
Standard Seating:	Compact molded plastic with cushion
Range Per Charge:	Up to 10 miles
Battery Charger:	Off-board, 2 amp
Battery Requirements:	(2) 12V deep cycle
Battery Size:	12 AH
Battery Weight:	9 lbs. each
Seat Weight:	17 lbs.
Heaviest Piece:	66 lbs.
Color:	Blue (only)

