


QUICKIE®

Q7

Beats Titanium...
Lighter Weight,
Lighter Price




QUICKIE Q7



Nothing Beats a Quickie™

The Quickie Q7 features 7000 Series Aluminum and ShapeLoc Technology that together shift the axis of the mobility world. Q7 is lighter than titanium, offers better performance and sleek style at an aluminum price.



Unparalleled Performance

- Weight Distribution
 - New caster wheel placement options make the chair easier to turn
- Rigidity
 - Ovalized tubing in axle plate and caster housing decreases frame flex
 - Increased frame bend strength results in better propulsion efficiency

Unique Style

- Hydroformed Caster Housing Design
- Color Options
 - New! Standard and Kolorfusion Frame Colors: Copper, Fuschia, Glow in the Dark, Rootbeer, Silver, Titanium, Zebra and Camo Patterns
 - Backrest Frame and Side Guards can be ordered to match



* Featured with 4-Link rear suspension

Unbelievably Lightweight

- Fully loaded configuration as low as 14 lbs.
- 7000 Series Aluminum
 - Lightweight Aluminum that can be tempered for strength through ShapeLoc technology



ShapeLoc Technology is the process of converting an aluminum frame into stronger, more rigid T6 Aluminum. The open frame is fabricated using 7000 Series Aluminum then heat treated to convert the metal to T6 Aluminum.

The Advantage: Stronger metal = thinner tube walls = an unbelievably lightweight chair!

SPECIFICATIONS Quickie Q7

Overall Weight	16" x 16": 17.5 lbs.
Weight Capacity	265 lbs.
Product Width	12" - 20"
Product Length	12" - 20"
Center of Gravity	1" - 4"
Back Height Adjustments	8" - 20"
Front Seat Height	16" - 21"
Rear Seat Height	13" - 20"
Frame Angle	75°, 80°, 85°
Warranty	Limited/Lifetime

QUICKIE Q7 www.QuickieQ7.com