

Hoyer® Pool Lift


- Great support for easier recreation and therapeutic swimming
- Stainless steel painted finish resists corrosion in humid environments
- Heavy-duty steel construction supports up to 400 lbs
- A Full 62" lifting range and 360° swivel helps gently place any client into water
- Hydraulic power can be used anytime



Specification	Imperial	Metric
Safe Working Load	400 lbs	181.8 kg
Length of Mast.....	57.5"	146 cm
Height of Mast Overall		
Cement-in Sleeve.....	51"	129.5 cm
Surface Mounted Base.....	59.25"	150.5 cm
Length of Boom	45"	114.3 cm
Height of Boom Minimum		
Cement-in Sleeve.....	24.25"	61.6 cm
Surface Mounted Base.....	32.25"	81.9 cm
Height of Boom Maximum		
Cement-in Sleeve.....	86.25"	219 cm
Surface Mounted Base.....	94.25"	239.4 cm
Lifting Range	62"	157.5 cm
Weight	59 lbs	26.8 kg
Packaged Weight	68 lbs	30.9 kg

Hoyer Sling Information

Model	Part Numbers	Material	Chains	Straps	Head Support	Commode	Warranty	HCPCS Code
112	112-D	Dacron	Yes	NA	No	No	6 months	E0621



112-D