

Ti

L I T E®



2GX Shown with Available Options



THE NEW TILITE 2GX.
TECHNOLOGY ENGINEERED TO EMPOWER.



THE NEW TILITE 2GX. DESIGNED TO INSPIRE!

The new 2GX titanium folding frame utilizes state of the art materials and design elements to maximize performance, function and durability. By combining TiLite's Tru-Fit component system and signature "bullet" caster housing with a unique curved axle plate, the 2GX is as adjustable as you need it to be. Test drive the new 2GX and experience the super smooth ride and performance advantage of the TiLite 2GX.

TILITE INNOVATION.

- "Tru-Fit" system.
- Rigid Folding Mechanism - "RFM".
- Adjustable "Bullet" caster housing.

TITANIUM. Ti 3-2.5 seamless aerospace grade titanium provides superior vibration dampening, reducing user fatigue and creating the 2GX's super smooth ride. Titanium also reduces fatigue on the wheelchair, resulting in greater durability and fewer repairs, in addition to providing superior comfort for the user.

SIT DOWN AND STAND OUT. The 2GX features TiLite's new Tru-Fit clamp-on component system; comprised of the RFM, rear axle plate and front caster clamp. These components all clamp onto the new 2GX frame instead of being bolted through. The Tru-Fit system results in a 50% reduction of holes in the frame, allowing for finite adjustments since there is no need to line up components with pre-drilled holes. This design innovation also greatly improves the durability of the 2GX frame.

RFM. TiLite's exclusive Rigid Folding Mechanism RFM cross tube design, first developed for the TX, allows TiLite to incorporate one-piece side frames into its newest line of folding chairs. This innovation increases performance by reducing flex because the side frames are welded together, not bolted. TiLite's RFM cross tubes also feature snap in seat saddles that lock into place, helping the 2GX achieve near rigid chair performance.

TECHNICAL SPECIFICATIONS

Frame Type	Titanium Folding
Frame Style	Fixed Front End - Full Adjustability Swing In or Out Removable Front End - Full Adjustability
Weight Limit	250 lbs.
Seat Width	12" - 20"
Seat Depth	12" - 20"
Front Seat Height	15" - 22"
Rear Seat Height	13" - 20.5"
Camber	0°, 2°, 4°, 6° or 8°
Front Angle	Fixed Front - 70°, 80° or 85° Swing Away - 70°, 80°, 85° or 90°
Front Wheels	18 Options
Rear Wheels	50 Options
Handrims	9 Options
Center of Gravity	1"- 6"
Chair Weight *	Fixed Front End - 13.5 lbs. Swing In or Out Removable Front End - 15.9 lbs.



PERFORMANCE



* Weights based on 16" x 16" wheelchairs with all standard features except rear wheels and wheel locks. Rear wheels are not included due to the varying weights of the numerous wheel configurations offered by TiLite. For reference, a pair of standard 24" rear wheels weighs 7.6 lbs. A pair of standard wheel locks weighs 0.9 lbs.



HAND MADE



STYLE