



THE ROHO[®] HARMONY[®] CUSHION

LIGHTWEIGHT

The HARMONY[®] Cushion from The ROHO Group provides your clients with the lightweight, stable solution they deserve. By combining a lightweight and comfortable ROHO AIR FLOTATION cushion with a pre-contoured foam base, encased within a fluid-resistant cover, your client can maximize their abilities, independence and freedom without the worry of adding unnecessary weight to their mobility system.

STABLE

REDUCES DISCOMFORT

Finding a cushion that reduces discomfort, matches your client's shape and addresses today's challenging funding climate shouldn't be difficult. Think HARMONY from The ROHO Group.

FUNDING- FRIENDLY



HARMONY[®] CUSHION

LIGHTWEIGHT

Don't limit client function with a heavy wheelchair cushion. Consumers will spend thousands of dollars to shave a few pounds off of a wheelchair, and then disregard how the cushion impacts the weight of the entire mobility system. A heavy mobility system can limit the client's function and increase the potential for long-term shoulder damage. At a lean 1.5 lbs. (1 kg), the HARMONY Cushion adds minimal weight to your client's mobility solution.




STABILITY

Stability builds a foundation of function. If a client is working to achieve balance or stability they may have difficulty achieving functional goals. The HARMONY cushion is built on an anatomically-shaped foam platform that enhances a client's seated stability and helps maximize their functional potential.

COMFORT

Many times comfort is dismissed as a luxury, but comfort and function go hand in hand. The HARMONY includes a ROHO AIR FLOTATION cushion positioned over the foam base that provides a shape-fitting interface reducing discomfort through enhanced blood flow and achieving longer sitting times.



APPLICATIONS	FEATURES & BENEFITS	MORE PRODUCTS FROM ROHO
<p>THE HARMONY CUSHION IS MOST EFFECTIVE FOR CLIENTS WHO:</p> <ul style="list-style-type: none"> - need a surface that reduces discomfort - require postural support - are at moderate to low risk for skin breakdown - are concerned about weight of mobility system 	<p>ROHO AIR FLOTATION CUSHION: Reduces Discomfort</p> <p>PRE CONTOURED FOAM BASE: Stability/Ease of Transfer</p> <p>LIGHTWEIGHT: Increased Mobility & Function</p> <p>WASHABLE: AIR FLOTATION cushion overlay may be cleaned with soap and water. Fluid-resistant machine-washable cover.</p>	 <p>CONTOUR SELECT[®] CUSHION</p>
<p>SPECIFICATIONS</p> <p>INCLUDES: HARMONY Cushion Base, AIR FLOTATION Cushion Overlay, cushion cover, operations manual, pump & patch kit.</p> <p>CONSTRUCTION: Specially formulated polyurethane die-cut foam base, polyvinylchloride AIR FLOTATION cushion overlay, flame and fluid-resistant cover.</p> <p>SIZES: Wide range to fit chair widths 14" to 20"/30cm to 51cm.</p> <p>WEIGHT: Approximately *1.5 lbs / 1 Kg. <i>*Varies by size.</i></p> <p>WEIGHT LIMIT: 250 lbs / 113 Kg.</p> <p>WARRANTY: 18-Month Limited Warranty.</p> <p>MEDICARE HCPCS BILLING CODE: E2607 Adjustable Combination Skin Protection and Positions Seat Cushion</p>	<p>We recommend consulting with a certified rehab professional when ordering and fitting ROHO Group products. Visit NRRTS.org for information on certified providers in your area.</p>	 <p>nexus SPIRIT[®] CUSHION</p>  <p>JETSTREAM PRO[®] BACK SUPPORT SYSTEM</p> <p>FOR ADDITIONAL INFORMATION, PLEASE CONTACT:</p> <div style="border: 1px dashed black; height: 100px; width: 100%;"></div>



100 North Florida Avenue • Belleville, IL 62221-5429 • USA • www.therohogroup.com
 U.S. : 1 800 851 3449 • 1 618 277 9173 • fax 1 888 551 3449
 Outside the U.S. : 1 618 277 9150 • fax 1 618 277 6518

©2009 ROHO, Inc. The following are registered trademarks of ROHO, Inc: ROHO[®], shape fitting technology[®], HARMONY[®], CONTOUR SELECT[®], JETSTREAM PRO[®], and nexus SPIRIT[®].

ROHO, Inc. has a policy of continual product improvement and reserves the right to amend specifications presented in this sell sheet.

LITG01020
9/09