



## THE ROHO<sup>®</sup> MOSAIC<sup>®</sup> CUSHION

### COMFORTABLE

For over 30 years, through research and design, The ROHO Group has manufactured products to fit the needs of a variety of users. The MOSAIC cushion is especially suited for users who require a cushion for basic level skin protection, support and comfort. Using air flotation technology, the MOSAIC cushion is not only comfortable, but the use of interconnected air cells to distribute forces evenly across the cushion allows it to be used effectively for low risk skin protection.

### COST EFFECTIVE

A cost effective cushion, the MOSAIC is designed to provide a comfortable, functional and stable sitting environment. MOSAIC is the affordable choice.

### LIGHTWEIGHT



# MOSAIC CUSHION

## COMFORTABLE

Flexible, interconnected air cells can be adjusted to closely shape the body's unique contours, providing customized comfort.



## COST EFFECTIVE

Constructed of polyvinyl, the MOSAIC cushion is an affordable solution for low risk skin protection needs in various facilities or home environment.

## LIGHTWEIGHT

The MOSAIC cushion weighs only 1 pound (.5 kg). Simple to transfer from one location to another; simple to use!



APPLICATIONS	FEATURES & BENEFITS	MORE PRODUCTS FROM ROHO
<p><b>THE MOSAIC CUSHION IS EFFECTIVE FOR CLIENTS:</b></p> <ul style="list-style-type: none"><li>- who require a basic level of skin protection</li><li>- who require comfort and support</li><li>- who have partial or full sensation in the seated area</li><li>- who have fair to good seated balance</li></ul>	<p><b>AIR FLOTATION TECHNOLOGY:</b> Provides comfort and is effective in skin protection.</p> <p><b>LIGHTWEIGHT:</b> Helps reduce fatigue by adding very little weight to the wheelchair or seating system.</p> <p><b>COST EFFECTIVE:</b> Constructed of polyvinyl, the MOSAIC cushion is an affordable solution for skin protection needs in various facilities or home environment.</p> <p>The MOSAIC is designed for clients who have partial or full sensation in the seated surface area. For clients who are at high risk or have existing skin breakdown, we recommend the use of ROHO DRY FLOTATION® Cushions.</p>	 <p><b>PRODIGY MATTRESS OVERLAY®</b></p>  <p><b>QUADTRO SELECT® HIGH PROFILE® CUSHION</b></p>
<p><b>SPECIFICATIONS</b></p> <p><b>INCLUDES:</b> MOSAIC cushion, cushion cover*, hand inflation pump, operating instructions and patch kit</p> <p><b>CONSTRUCTION:</b> Polyvinyl Chloride</p> <p><b>SIZES:</b> Available to fit chairs: 16" x 16"/41cm x 41cm 18" x 16"/46cm x 41cm 18" x 18"/46cm x 46cm</p> <p><b>WEIGHT:</b> Varies by size, approximately 1 lb./ .5kg (based on 18" x 16" size)</p> <p><b>WEIGHT LIMIT:</b> 250 lbs. / 114 kg</p> <p><b>WARRANTY:</b> 12-month limited</p> <p><b>MEDICARE PART B CODE:</b> E2601 General Use Wheelchir Seat Cushion</p> <p>*The MOSAIC cushion is available with a standard cover or an incontinent cover. Please specify which model you prefer upon placing your order.</p>		<p><b>FOR ADDITIONAL INFORMATION, PLEASE CONTACT:</b></p> <div></div>



100 North Florida Avenue • Belleville, IL 62221-5429 • USA • [www.therohogroup.com](http://www.therohogroup.com)  
U.S. & Canada: 1 800 851 3449 • 1 618 277 9173 • fax 1 888 551 3449  
Outside the U.S. & Canada: 1 618 277 9150 • fax 1 618 277 6518

©2007, 2009 ROHO, Inc. The following are registered trademarks of The ROHO Group: ROHO®, DRY FLOTATION®, MOSAIC®, PRODIGY Mattress Overlay®, and QUADTRO SELECT®.

ROHO, Inc. has a policy of continual product improvement and reserves the right to amend specifications presented in this sell sheet.

LIT602099  
09/09