

Invacare® Tracer® IV wheelchair features superior durability, rollability and streamlined looks. It is designed for the active individual who needs the strength of a reinforced heavy-duty frame. Its many options include components interchangeable within the Tracer series. The Tracer IV wheelchair is long-lasting and low-maintenance.

Features

- Dual axle position allows variation of seat-to-floor heights
- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Reinforced frame includes gusseted sides, caster journals and dual crossbraces
- Corded urethane rear tires (350 lb.) and pneumatic flat free tires (450 lb.) mounted on “no flex” wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Upholstery is dual-embossed



Invacare® Tracer® IV Wheelchair



Heavy-duty package increases weight capacity to 450 lb. and consists of heavy-duty composite wheels with pneumatic tires with flat-free inserts and semi-pneumatic casters.

Specifications

Seat sizes		Arm styles	
Widths	18", 20", 22", 24"	Full-length fixed height, padded, removable	
Depths	18" (standard), 20" (by seat extension)	Desk-length fixed height, padded, removable	
Seat-to-floor heights	17.5", 19.5"	Full-length adjustable height, padded, removable	
Back height	16" fixed	Desk-length adjustable height, padded, removable (All arms are conventional)	
Casters and wheels		Weight capacity	350 lb., 450 lb.*
Casters	8"	Product weight (T4/WD28/28)	42 lb.
Wheels	24"	Limited warranty	5 years on frame and crossbraces
Front riggings		* Heavy-duty package only.	
Swingaway footrest	17" - 21" knee-to-heel		
Swingaway padded elevating legrest	17.5" - 21" knee-to-heel		
Articulating elevating legrest	17" - 25" knee-to-heel		
(Footplates available in aluminum and composite)			