



**Invacare® Tracer® SX5 wheelchair's** lightweight frame, weighing less than 36 pounds, makes this chair perfect for rental and long- or short-term use. This wheelchair offers multiple seat widths, arm styles, sectional back and seat-to-floor heights by utilizing the dual axle position frame base. The 20" and 22" widths include a semi-pneumatic caster to accommodate a weight capacity of 300 lb. Combining the design and technology of the Invacare 9000 series wheelchairs makes it possible to interchange components with the Tracer SX5 wheelchair.

Invacare® Tracer® SX5 Chair

**Features**

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on “no flex” wheels offer superior performance
- Sectional back available in 16" and 18" back heights
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Meets RESNA standards for manual wheelchairs in product performance

**Specifications**

Seat sizes	
Widths	14", 16", 18", 20", 22"
Depths	16"
Seat-to-floor heights	
	17.5", 19.5"
Back heights	
	16", 18" sectional
Casters and wheels	
Casters	6", 8"
Wheels	20", 22", 24"
Front Riggings	
Hemi swingaway footrest	17"-21" knee-to-heel
Hemi swingaway padded elevating legrest	17.5"-21" knee-to-heel
Articulating elevating legrest	17"-25" knee-to-heel

(Footplates available in aluminum or composite)

Arm styles	
	Full-length fixed height, padded
	Desk-length fixed height, padded
	Full-length adjustable height, padded
	Desk-length adjustable height, padded
(All arms are flip back non-removable)	
Weight capacity	
	250 lb. (300 lb. on 20" and 22" width)
Product weight ( TRSX5/WD86/28FB)	
	34 lb.
Limited warranty	
	5 years on frame and crossbraces