

# KN-800B

## Standard fixed-arm wheelchairs



- Folding powder coating steel frame
- Padded flip back armrests
- Swing-away detachable footrests
- 24" x 1" mag wheels with rubber tires
- 8" x 1" caster wheels with rubber tires
- Adjustable toggle wheel locks.
- Adjustable length footplates.
- Side panels and anti-slip footplates
- High-grade flame retardant embossed

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KN-800B	18	16	20	16	36.5	25	37	250
KN-800B-N	16	16	20	16	36.5	23	37	250

Model No.	Description
KN-800B	18", fixed arms, swing-away footrests
KN-800B-E	18", fixed arms, elevating legrests
KN-800B-N	16", fixed arms, swing-away footrests
KN-800B-N-E	16", fixed arms, elevating legrests

Common Options: \*Not available on all models \*\* Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Quick Release Axle • Spoke Wheels • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335