

GetGo Classic Rolling Walkers

NOVA



4202CBL - Blue



4202CRD - Red



4202CGN- Green

**Be fashionably pink
with the new 4202CPK!**



Specifications

Weight Capacity:	300lbs
Weight:	13.5 lbs
Handle Ht Adj:	32.5" - 36.5"
Seat height:	22"
Seat Dimensions:	13.5"w x 12.5"d
Width btw Handles:	17.5"
Overall Dimensions:	24"w x 23"d
Approx Patient Ht:	5'4"-6'1"
Overall Dimensions	
Folded:	36"l x 24"w x 9.5"d
Bag Dimensions:	15"w x 9"l x 6" d
Wheel Size:	6"

Nova Ortho-Med, Inc.
1470 Beachey Pl
Carson, CA 90746



**For Customer Service, call
Toll-Free at (800) 557-6682**