

BLUE STREAK



Blue Streak Single Axle (K0001)

FEATURES:

- Blue powder coated frame, black cross brace.
- Flip back, padded desk length arm.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Solid rubber tires mounted on composite wheel provides durability and low maintenance.
- Comes with push to lock brakes.
- Comes with calf strap.
- Swing away footrests with plastic footplates come standard.
- Color: Black seat and back, blue frame.
- 3 year limited warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*