

VIPER PLUS FULL RECLINING



**NEW STATE-OF-THE-ART
HYDRAULIC RECLINING
MECHANISM.**

Viper Plus Reclining Deluxe High Strength, Lightweight (K0001/K0007)

FEATURES:

- New state-of-the-art hydraulic reclining mechanism.
- Carbon steel frame with double coated chrome for an attractive, chip-proof, maintainable finish.
- Durable, heavy gauge flame retardant nylon upholstery resists mildew and bacteria.
- Nylon, padded armrests provide added patient comfort.
- Flip-back arms allow for easy transfer.
- Composite Mag-style wheels are lightweight and maintenance free.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Standard dual axle provides easy transition of seat to hemi-level.
- Composite footplates are durable and lightweight.
- Front caster forks are adjustable in 3 positions.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- Comes with push to lock wheel locks.
- Comes standard with swing-away elevating legrests.
- Color: Black
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*