

CIRRUS PLUS EC



**COMES
STANDARD WITH
SEAT BELT AND
SEAT CUSHION**



**BACK POCKET
FOR EASY AND
CONVENIENT
STORAGE**

Cirrus Plus EC Folding Power Chair

CPN16FBA	16" Wheelchair, with Flip Back Arms, 1/cs
CPN18FBA	18" Wheelchair, with Flip Back Arms, 1/cs
CPN20FBA	20" Wheelchair, with Flip Back Arms, 1/cs
CPN16FBA-POM	16" Wheelchair, with Flip Back Arms, (Comes with "Peace of Mind Service"), 1/cs
CPN18FBA-POM	18" Wheelchair, with Flip Back Arms, (Comes with "Peace of Mind Service"), 1/cs
CPN20FBA-POM	20" Wheelchair, with Flip Back Arms, (Comes with "Peace of Mind Service"), 1/cs

FEATURES:

- Foldable, lightweight frame makes portability and transportation easy and convenient.
- Integrated, programmable PG controller.
- Adjustable height flip back and detachable arms.
- Attractive carbon steel frame with silver vein finish.
- Flat free, 12.5" black drive wheels.
- Front 8" flat free black casters.
- Anti-tippers standard.
- Adjustable tension back upholstery.
- Comes standard with seat belt, seat cushion with pouch, swing-away footrest with calf straps and heel loops.
- Elevating legrests optional. (Item # LELR-TF)
- Folding battery tray.
- Padded back upholstery and cushioned seat.
- 22" model has double cross brace frame.
- Adjustable length joystick mount.
- Wheel locks.
- Comes with batteries.
- Has been tested and surpassed the ANSI RESNA testing requirement.
- Attractive sling style seating.
- 300 lb. Weight Capacity.
- 400 lb. weight capacity.



**PWC
GROUP 2 HD
HCPCS
CODE:
K0824**

Cirrus Plus HD Folding Power Chair

CPN22FBA	22" Wheelchair, with Adjustable Height/Width Removeable Arms, 1/cs
CPN22FBA-POM	22" Wheelchair, with Adjustable Height/Width Removeable Arms, (Comes with "Peace of Mind Service"), 1/cs

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.