



TR17



CONVENIENT CUP HOLDER
CAN EASILY ATTACH TO CHAIR!
SOLD SEPARATELY MODEL # STDS1040S

DURABLE
CHROME FINISH!



COMES WITH
TOOL FREE ADJUSTABLE
FOOTRESTS

Deluxe Go-Kart Steel Transport Chair (Chrome)

TR17	17" Chrome Frame, Black Upholstery, 1/cs
TR19	19" Chrome Frame, Black Upholstery, 1/cs
LELR-TF	Swing-away Elevated Legrest, Tool Free, 1pr/cs, (Optional front riggings)

FEATURES:

- Weighs under 26 lbs.
- Available in 17" and 19" seat widths.
- Carbon steel frame with double coated chrome provides an attractive, chip-proof, easy to maintain finish.
- Aluminum casted back release hinge allows back to fold down for easy storage and transport.
- Comes with carry pocket on back rest.
- Back folds down for easy storage and transport.
- Comes with tool free swing-away footrests.
- Comes with seat belt for added safety.
- 8" casters in front and rear.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Rigging	Product Weight	Carton Shipping Weight	Weight Limit
TR17	20"	9"	36"	19"	18"	17"	16"	33"	26 lbs.	31.5 lbs.	300 lbs.
TR19	22"	9"	37.25"	19"	18"	19"	16"	33"	26 lbs.	31.5 lbs.	300 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

transport chairs

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM