

# BARIATRIC SENTRA EXTRA HEAVY DUTY

COMES WITH  
DUAL CROSS  
BAR!

AVAILABLE  
IN 20", 22"  
AND 24"  
SEAT  
WIDTHS.

NOW  
COMES WITH  
RIGGINGS IN  
CARTON!



SHOWN WITH OPTIONAL STDELRTF  
SWING-AWAY ELEVATING LEGRESTS.

## Sentra Heavy-Duty

Dual Axle when used with STDS1074 Caster  
(Optional Caster Sold Separately) (K0007)

### FEATURES:

- Carbon steel frame with TRIPLE coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable, heavy gauge, reinforced naugahyde upholstery.
- Padded armrests provide added patient comfort.
- **Composite Mag-style wheels with chrome handrim are lightweight and maintenance free.**
- Dual position front fork.
- Dual position front frame.
- Comes standard with carry pocket on back rest.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Aluminum footplates are attractive, crack-proof and lightweight.
- **"Hemi" adjustable caster is optional (Part # STDS1074).**
- **Reinforced side frame, caster journals and dual cross bars.**
- Comes with push to lock wheel locks.
- Black upholstery, chrome frame.
- Limited Lifetime Warranty.

BARIATRIC

**500 LB**

WEIGHT CAPACITY

*Due to product improvements and changes, all images in this catalog are subject to change.*

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.