



Roscoe Mobility Transport Chair

The Roscoe Transport Chair has a carbon steel frame with a double-coated chrome or powder-coated finish that is attractive, chip-proof, and easy to maintain. Weighing less than 26 pounds, it is equipped with comfortable, padded nylon and a folding backrest that makes it easy to pack and go. Comes complete with a seat belt for added safety. Tool-free swing-away, aluminum footrests with heel loops come standard.

Product Features

- Carbon steel frame with double-coated chrome or powder-coated finish that is attractive, chip-proof, and easy to maintain.
- Permanent rear axle
- Black padded nylon upholstery is lightweight and easy to clean
- Comes with upholstery handles to fold chair
- Back release hinge allows back to fold down for easy storage and transport
- Fixed, full-length, padded arms provide additional comfort
- Comes with seat belt for added safety
- Swing-away, aluminum footrests with heel loops come standard
- Composite 8" x 1" wheels with solid rubber tires are lightweight and maintenance free
- Easy push-to-lock rear wheel locks
- Limited Lifetime Warranty

Product Specifications

- Product Weight: 25 lbs.
- Overall Width with Arm: 23½"
- Arm Style: Fixed
- Overall Depth: 27" (without front rigging)
- Seat Width: 19"
- Seat Depth: 16"
- Seat-to-Floor Height¹: 19"
- Back Height: 16"
- Back Style: Folding
- Weight Capacity: 250 lbs.
- HCPCS Code: E1038*

Item #	Finish	Seat Width	Armrests	Front Riggings
KT1916FLFSA	Chrome	19"	Fixed Full-Length Arms	Swing-Away Footrests
KT19BG	Burgundy	19"	Fixed Full-Length Arms	Swing-Away Footrests
KT19BL	Blue	19"	Fixed Full-Length Arms	Swing-Away Footrests

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice. ¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.



To Order:
(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com