

SilvarBoom 250

Base weight 58 lbs. with Straight Arm¹ 65 lbs. with Offset Arm¹ Maximum lifting capacity 250 lbs. Exterior rated Yes Horizontal strap adjustment 21" - 35"2 Boom arm angleadjustability 0°, 5°, 10°, 12°, 15° Motor 12 volt DC Warranty 3-year transferable Vehicle Platform Full-size van. Minivan, Truck, Sedan, SUV and Crossover

LIFTS CURRENT PRIDE® PRODUCTS³

All Pride® Scooters4, except: Pursuit® XL & Wrangler®

All Pride® & Quantum® PWCs⁵, except:
Q6 Edge®, Q6000Z, Q6000Z HD, Q640, R-4000, R-4400, Q1121, Q1420, Q1424, Q1650, Jazzy® 1170 XL Plus, Jazzy® 1420 & Jazzy® 1650




Call for details or visit "Find My Compatible Lift" at prideliftsandramps.com for more applications.

- Includes post weight
- Adjustable in 1" increments
- Vehicle dependent
- Pursuit® can be lifted with seat off only
- 5 Seating system dependent

Note: Please contact your Pride® Representative for compatibility of other manufacturers' mobility devices.

SilverBoom 250







www.prideliftsandramps.com • www.pridemobility.com (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)





www.prideliftsandramps.com • www.pridemobility.com (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)



www.prideliftsandramps.com • www.pridemobility.com (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

SIIvarBoom 250

The SilverBoom 250 easily fits into most vehicles and is a perfect complement to a Jazzy® Power Chair or Pride® Scooter. Engineered to provide smooth, efficient operation, maximum durability, and a sporty style, the SilverBoom 250 makes it easy and affordable to take your mobility product with you wherever you go.

- Exterior rated
- Easy-to-use hand control
- Adjustable to accommodate a variety of applications
- Controller with on/off switch
- Straight or Offset Arm available
- Compatible with dual-charge battery box
- Lifting capacity 250 lbs.
- · Designed to prevent strap twisting
- · Stop warning on strap to identify lifting limit
- Offers angle adjustability at 0°, 5°, 10°, 12° and 15°
- Boom arm allows horizontal strap adjustment to reach from 21" to 35" from base
- Lift ships fully assembled
- Custom drill and no-drill installation bases available

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Interior/Exterior Lift

Straight Arm Application

















Pride Lifts & Ramps