

Backpacker[®] AVP

Base weight lift 193 lbs.¹

Scooter/power chair overall weight 350 lbs.

Maximum product width 27"

Maximum travel distance from the ground 37"

Motor Sealed 12 volt DC

Overall platform length 42" - 34"

Maximum product wheelbase 41"

Warranty Three-year transferable

Optional battery pack 12 volt, on-board

Optional charger 110 volt

¹ Weight does not include any accessories

² Vehicle dependent

Can easily lift the following current Pride[®] Power Chairs²

- Jazzy[®] 600
- Jazzy[®] 614/614 HD
- Jazzy[®] 1113 ATS
- Jazzy Select[®]
- Jazzy Select[®] 6
- Jazzy Select[®] 6 Ultra
- Jazzy Select[®] Traveller
- Jazzy Select[®] GT
- Jet[™] 3 Ultra
- Go-Chair
- Z-Chair

Can easily lift the following current Quantum[®] Power Chairs²

- Quantum[®] 600/600 XL/640
- Quantum[®] 610
- Quantum[®] 614
- Quantum[®] 1103 Ultra
- Quantum[®] 1121
- Quantum[®] 6400
- Quantum[®] 6000Z/6400Z/6000Z HD
- Quantum[®] R4000/R4400

Can easily lift the following current Pride[®] Scooters²

- Celebrity[®] X 3 & 4
- Go-Go[®] Ultra X 3 & 4
- Go-Go Elite Traveller[®] Plus 3 & 4
- Go-Go Elite Traveller[®] Plus HD 3 & 4
- Legend[®] 3 & 4
- Maxima 3 & 4
- Revo[®] 3 & 4
- Victory[®] 9 3 & 4
- Victory[®] 9 PS
- Victory[®] 10 3 & 4

¹ Weight does not include any accessories

² Vehicle dependent

Battery Pack

NEW, OPTIONAL MULTI-PURPOSE BATTERY BOX



- Industry first
- Industry exclusive
- Charges the mobility device from the lift as you drive

Backpacker[®] AVP



Silver Star[®] Lifts & Ramps

www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

Backpacker® AVP

Now you can easily take your mobility on the road in your sport utility vehicle with the Backpacker® AVP from Silver Star® Lifts and Ramps. It's the ultimate solution to lifting your powered mobility device into your vehicle; engineered specifically to be a perfect fit for a wide variety of popular vehicles including SUVs. Like all Silver Star Lifts, the Backpacker AVP features all-steel construction and a protective finish, along with an industrial rated 12 volt motor and battery pack, and easy-to-use controls to combine rugged durability with safe, simple and reliable operation.

Note: The Backpacker® AVP can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting with a power positioning system. It can fit into most SUVs and mini-vans with little or no vehicle modification; full-size vans require modification.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



- Adjustable platform for smaller SUVs 42"-34"
- Removable, front drive panel for easier serviceability
- New, adjustable wheel chocks on the platform surface
- Single piece horizontal drive rail: minimum welding points for noise reduction, increased rigidity and more fluid movement
- New, angled low-profile platform
- Easy-to-use two button hand control
- Can lift a wide range of Pride® Scooters and Power Chairs and other manufacturers' mobility devices

- 350 lbs. lifting capacity
- Three-year transferable warranty
- Fits into SUVs with a 31" opening height
- Manual override
- Improved micro switch protection
- False floor available for specified vehicles through special order to preserve vehicle aesthetics
- Ratcheting straps now available as an option