



**Pride® Lift Chairs**  
*Experience the Luxury*

Comfort and style complement one another perfectly in Pride's Heritage Collection Lift Chairs. With an array of plush fabrics to fit any décor and a superior design that provides exceptional strength, Heritage Collection chairs deliver style and durability that will last for years to come. Plus, Pride's Heritage Collection offers numerous standard features, superb comfort, and quiet and smooth performance, all at a great value.

# HERITAGE COLLECTION

## HERITAGE COLLECTION MODELS



GL-58	GL-358PW	GL-358S	GL-358M	GL-358L	GL-358XL	GL-358XXL
2-Position Medium 5'2" - 5'10"	3-Position Small Wide Chaise Lounger 5'3" - Below	3-Position Small Chaise Lounger 5'3" - Below	3-Position Medium Chaise Lounger 5'4" - 6'0"	3-Position Large Chaise Lounger 5'9" - 6'2"	3-Position Heavy-duty Chaise Lounger 5'6" - 6'1"	2-Position Very Heavy-duty Chaise Lounger 5'6" - 6'1"
					<b>500 LBS. WEIGHT CAPACITY</b>	<b>600 LBS. WEIGHT CAPACITY</b>

## SPECIFICATIONS


	GL-58	GL-358PW	GL-358S	GL-358M	GL-358L	GL-358XL	GL-358XXL
<b>Seat To Floor:</b>	20"	19"	19"	20"	21"	20"	20"
<b>Seat Depth:</b>	20"	19"	18"	20"	22"	21"	22"
<b>Seat Width:</b>	20"	22"	18"	20"	22"	26"	30"
<b>Top of Back to Seat:</b>	24.5"	23.5"	23.5"	24.5"	25.5"	29"	29"
<b>Distance From Wall:</b>	11.5"	17"	17"	17"	17.5"	18"	18"
<b>Overall Chair Width:</b>	35"	39"	35"	35"	37"	39.5"	43"
<b>Chair Weight:</b>	108 lbs.	118 lbs.	118 lbs.	117 lbs.	120 lbs.	159 lbs.	168 lbs.

## FABRIC SELECTION


### NEW STANDARD FABRICS



### OPTIONAL FABRICS



### PADDED-SUEDE



### MICRO-SUEDE



### LEVANTE VINYL




### NEW CRYPTON® SUPER FABRICS - SMART SUEDE



### NEW PLUSH VELVET COLLECTION



### NEW WINDSOR MICRO-SUEDE COLLECTION



### VINYL



## BENEFITS

- Engineered furniture grade laminate/hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Standard head and arm covers
- Fire retardant (California specifications; bulletin 117)
- Standard, large dual pockets
- Removable back - no tools required
- 375 lbs. weight capacity (GL-358XL supports up to 500 lbs.; GL-358XXL supports up to 600 lbs.)
- Integrated battery backup in case of power failure<sup>1</sup>
- Transformer features self-diagnostic electronics with 12' cord
- Lifetime warranty on steel lift frame
- 7-year pro-rated warranty

## OPTIONS\*

- Heat & Massage
- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Leather/vinyl wipes

Additional lead-time may apply for some options.

\* Additional charges apply for all options.

<sup>1</sup>. Batteries not included.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.