



Roscoe Mobility Rollator with Padded Seat

The Roscoe Rollator with soft grip tires is perfect for indoor and outdoor use. The rollator's handles adjust to the user's height while the reliable, loop style brakes control speed and lock wheels when seated. Comes equipped with a handy wire basket and a straight padded backrest. The flip-up padded seat features a zipper compartment for storage and an easy folding mechanism.

Product Features

- Product Weight: 16 lbs.
- 6" Wheels with Wheel Locks
- Loop Style Brakes
- Padded Seat with Easy Folding System
- Straight Padded Backrest
- Handy Wire Basket for Storage
- Limited Lifetime Warranty

Product Specifications

- Seat Measurements: 14½"(w) x 13½"(d)
- Seat-to-Floor Height: 23 ¼"
- Adjustable Handle Height: 32 ½" - 37"
- Unit Width: 20"
- Weight Capacity: 275 lbs.
- HCPCS Code: E0143*
- HCPCS Code - Seat: E0156*

Item #	Description	Color	Case Qty
ROS-RL10040A-BL	Rollator with Padded Seat, 6" Wheels	Blue	1
ROS-RL10040A-RD	Rollator with Padded Seat, 6" Wheels	Burgundy	1
ROS-RL6GR	Rollator with Padded Seat, 6" Wheels	Granite Gray	1

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com