

Pulse Oximetry

Six-Display Mode Pediatric Pulse Oximeter

The Roscoe Medical Six-Display Mode Pediatric Pulse Oximeter features a child-friendly design to help eliminate anxiety while monitoring oxygen level and heart rate. Its compact design includes a special mark to help children with proper finger placement. Low power consumption allows the unit to run for 30 hours of continuous operation with two AAA batteries. The all-in-one integrated device has six display modes and ten levels of adjustable brightness. One year warranty.



Attractive Retail Packaging!



Product Features

- Child-friendly design helps to reduce anxiety
- Special mark on unit helps children with proper finger placement
- Simple operation for immediate real-time spot checks
- Low power consumption allows for more than 30 hours of continuous operation with two AAA batteries
- Large, easy-to-read, dual color OLED screen with adjustable ten levels of adjustable brightness
- Six display modes
- Unit automatically powers off to save battery life when not in use for 8 seconds
- Low voltage warning indicator light
- Includes two AAA batteries, operating manual and lanyard
- One year warranty

Product Specifications

- Display Type: OLED
- SPO2 Measurement Range: 70 - 99%
Accuracy: 80% - 100%, $\pm 2\%$; 70% - 79%, $\pm 3\%$; $\leq 69\%$ no definition
- Pulse Rate Display Range: 30 - 235 BMP
Accuracy: 30 ~ 99 bpm, ± 2 bpm; 100 - 235 bpm, $\pm 2\%$
- Pulse Rate Display Mode: bar graph Indicator
- Data Update Time: < 15 s
- HCPCS Code: E0445*

Item #	Description	Qty
50020	Six-Display Mode Pediatric Pulse Oximeter	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com