

PostureLITE™

BY ROHO®

Simple. Adaptable. All the Time.

Air and foam cushion insert combines to self-adjust to your comfort.

Fluid-resistant cover is easily removed for cleaning.

Cushion cover's nonskid bottom materials keeps PostureLITE in place.

Perimeter air expansion chamber self-adjusts to conform to the contours of the body.

Performance In A Simple Package

PostureLITE™ Cushion by ROHO® is a new self-adjusting wheelchair cushion that combines unparalleled ease of use with comfortable skin protection. It is ideal for anyone who may be at risk of skin breakdown and may have normal sensation.

The preset air and foam cushion offers a perimeter air expansion chamber that conforms to the contours of the body to help provide comfortable suspension and skin protection. Similar to ROHO cushions, PostureLITE works by allowing continuous redistribution of air within the cushion with every weight shift. This maximizes the contact area of the cushion to help manage pressure and increase tissue viability.

The super-lightweight cushion weighs only 1.5 lbs. (0.6 kg) making it one of the lightest cushions available on the market.

Features

PostureLITE features a “closed system” design that utilizes both foam and air to provide fail-safe support which never needs adjustment. The flexible perimeter air expansion chamber allows the PostureLITE to adapt to various chair sizes making transferring the cushion to other seating surfaces simple and easy. The cover is a stretchable, fluid-resistant material that is machine washable.

Construction

Foam insert inside air bladder: Polyurethane foam.

Preset, sealed air bladder with perimeter baffles: PVC (DEHP and BPA free).

Cover: Top and sides are Polyurethane-coated Polyester. Nonskid bottom is PVC (DEHP and BPA free).

Easy-Fit Sizes

Three sizes adapt to fit a wide range of wheelchairs, from 15 to 20 inch widths (38 cm to 51 cm):

POSTUREL1: 16.75 in x 17.25 in (42.5 cm x 44 cm)

POSTUREL2: 18.75 in x 17.25 in (47.5 cm x 44 cm)

POSTUREL3: 20.00 in x 17.25 in (51.0 cm x 44 cm)

Cushion Height

3 in (7.5 cm) without load.

Weight

Varies by size, approx. 1.5 lbs (0.6 kg) – based on POSTUREL2 size.

Weight Limit

350 lbs (158 kg).

Warranty

12-Month Limited Warranty.

Certifications

Complies with EN 1021-1, EN 1021-2 and CAL 117 flammability requirements.

Funding Qualifications

U.S. – E2601 Code Verified. Consult your country's ROHO distributor for funding information.



PostureLITE™ gives you:

COMFORT

Self-adjusting design provides individualized comfort

CONVENIENCE

Requires no setup or adjustments

PROTECTION

Air redistribution design assists in the process of skin protection

DURABILITY

High quality materials, superior design and construction

EASY FIT

Sizes compatible with most wheelchairs

ENJOY A
NEW LEVEL
LIGHTWEIGHT
PROTECTION!



*“You will LOVE
the Simplicity
and Comfort!”*



PostureLITE™
BY ROHO®

100 North Florida Avenue • Belleville, IL 62221-5429 • USA

www.therohogroup.com

U.S.

Telephone: 1 800 851 3449 • 1 618 277 9173

Fax: 1 888 551 3449

Outside the U.S.

Telephone: 1 618 277 9150

Fax: 1 618 277 6518

© 2011 ROHO, Inc. PostureLITE™ is a trademark of Comfort Concepts, LLC.
The following is a registered trademark of ROHO, Inc.: ROHO®.
PostureLITE is distributed by ROHO, Inc.
ROHO, Inc. has a policy of continual product improvement and reserves the right to amend specifications presented in this document. The current version of this document is available on our web site, www.therohogroup.com.

LITG02168 Rev 9/11
2M
SPC 9/11

Introducing!

PostureLITE™
BY ROHO®

Innovations in Seating Technology by **ROHO®**