



Jazzy[®]

ELITE 6

The Jazzy[®] Elite 6 provides stability combined with the advantages of in-line motor technology. Front and rear casters offer solid performance indoors and out.



Jazzy[®] Elite 6
shown with Jazzy[®]
Armor™ (standard)

Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	2"
Turning radius	25"
Overall length	36.5" without front riggings
Overall width	23.25"
Drive wheels	9"
Front wheels	5" solid casters
Rear wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 13 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	21 lbs. each
Base weight	73.5 lbs.
Standard seat weight	37.5 lbs. (comfort high-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy Elite 6-2S-SS	(Code: K0822)
Jazzy Elite 6-2S-C	(Code: K0823)

FEATURES

- Casters on front and rear for stability
- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Comfort high-back seat with headrest
- New Jazzy® Armor™ includes shroud and controller guards to protect against daily wear and tear.

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Pride Lifts & Ramps

- Backpacker® Plus
- Outlander
- Backpacker® MV
- PrideBoom 250



COLOR CHOICES



Red



Blue

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

(US) 800-800-8586 ▪ (Canada) 888-570-1113 ▪ www.pridemobility.com