

The AmeriGlide Traveler brings exceptional value and easy transporta tion together in a compact, lightweight package. The stylish Traveler has a weight capacity of 275 lbs. and disassembles quickly and easily.



RAVELER

Specifications	
Max. weight capacity	275 lbs.
Maximum speed ¹	Up to 4 mph
Ground clearance ²	2.5″ at mid-frame
Turning radius ²	32.75″
Overall length ²	38″
Overall width ²	19.5″
Tires	
Front	8" solid
Rear	8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 1.5 amp
Per-charge range ^{1,5}	Up to 6 miles
Battery requirements ³	12V deep-cycle
Compatible batteries	12 AH
Weight	
Batteries (2 required)	12 AH (16 lbs. battery pack) 4
Weight w/o battery pack ²	72.5 lbs.
Heaviest piece ²	27.5 lbs. (front section)
Warranty	5-year limited on frame; 2-year limited on electronics; 1-year limited on motors
Standard colors	

Standard colors

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of plus or minus 10%.

- ² Due to manufacturing tolerances and continual product improvement, this
- specification can be subject to a variance of plus or minus 3%.
- ³ AGM or Gel-Cell type required.
- ⁴ Battery weight may vary based on manufacturer.
- ⁵ This range was obtained with the 4-wheel model.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

Features

Easily disassembles into 5 super lightweight pieces for convenient transport and storage

Sleek sporty style

Modular design for easy serviceability

Front frame-mounted seat post offers maximum stability

PDAC coded K0800 power operated vehicle



Easy Disassemble

