

| | |
|---|---|
| Model | Go-Go® ES 2 ¹ |
| Part Number | S81 |
| Weight Capacity | 250 lbs. |
| Maximum Speed | Up to 4 mph |
| Ground Clearance ⁴ | 1.5" at motor / 0.75" at anti-tips |
| Turning Radius ⁴ | 31" |
| Overall Length ⁴ | 36.75" |
| Overall Width ⁴ | 19.5" |
| Seat-to-Ground Height ³ | 20"- 21" |
| Seat-to-Deck Height ³ | 16.25"-17.25" |
| Front Tires | 2" x 7.5" solid |
| Rear Tires | 3" x 7.5" solid |
| Range Per Charge ^{2,3} | Up to 6.24 miles |
| Total Weight | With batteries: 77.5 lbs. Without batteries: 58 lbs. |
| Heaviest Piece When Disassembled ⁴ | 41.4 lbs. (scooter base) |
| Standard Seating | Type: Foldable, compact, vinyl-covered Weight: 17 lbs. Material: Black vinyl Dimensions: Width: 17" Depth: 15" |
| Motor Output | 24-volt DC motor, sealed mini-transaxle |
| Battery Requirements ⁵ | Size: (Two) 12 AH Weight: 9 lbs. each |
| Battery Charger | Off-board, 2 amp |
| Color | Black |

HCPCS Code: K0800 (Group 1, Standard)

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁴ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁵ Battery weight may vary based on manufacturer.

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

0003693-110821

GO-GO® es 2 ¹

OUTSTANDING VALUE...

...EXCEPTIONAL DURABILITY



FDA Class II Medical Device¹

LIVE YOUR BEST®

pridemobility.com |    

OUTSTANDING VALUE...

...EXCEPTIONAL DURABILITY

Foldable, compact,
vinyl-covered
seating with
armrests

Modular design
for easy
serviceability



Two speed
settings for
indoor and
outdoor use

Removable
front canvas
basket

The Go-Go® ES 2 from Pride® combines exceptional value and easy transportation together in a compact, lightweight package. The stylish Go-Go ES 2 has a weight capacity of 250 lbs., variable motor speeds and a convenient canvas basket located on the tiller.

Additional Features

- Lightweight, one-piece frame for added stability and easy portability
- Modular design for easy serviceability

Color



Black

Accessories

- Weather cover
- Military patches



Large, removable front
canvas basket.

High and low speed settings for
indoor and outdoor use.

Foldable, compact, vinyl-covered
seating with armrests.