

# LC-525 SERIES



The LC-525 Lift Chair Series features superior comfort, performance and style. A “zero-gravity” position, infinite-position hand control and super soft fabric selections make the LC-525 the ultimate in lift chair comfort.

### Features and Benefits:

- Tilt-in-space “zero-gravity” position with active seat box for unmatched comfort
- One-button push to “zero-gravity” position
- True Trendelenburg position
- Available in Petite Wide, Small, Medium and Large




	<b>LC-525PW</b>	<b>LC-525S</b>	<b>LC-525M</b>	<b>LC-525L</b>
	Infinite-Position “Zero-gravity” Chaise Lounger 5’2” & Below	Infinite-Position “Zero-gravity” Chaise Lounger 5’3” & Below	Infinite-Position “Zero-gravity” Chaise Lounger 5’4” - 6’0”	Infinite-Position “Zero-gravity” Chaise Lounger 5’9” - 6’2”

<b>Weight Capacity</b>	375 lbs.	375 lbs.	375 lbs.	375 lbs.
<b>Seat-to-Floor</b>	19”	19”	20”	21”
<b>Seat Depth</b>	19”	18”	20”	22”
<b>Seat Width</b>	24”	18”	20”	22”
<b>Top of Back to Seat</b>	25.5”	25.5”	26.5”	27.5”
<b>Back Style</b>	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
<b>Fully Padded Chaise</b>	Yes	Yes	Yes	Yes
<b>Distance From Wall</b>	28”	28”	29”	30”
<b>Overall Chair Width</b>	38”	32”	33”	35”
<b>Chair Weight</b>	130 lbs.	123 lbs.	126 lbs.	128 lbs.
<b>Warranty</b>	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated

### AVAILABLE FABRICS

#### CLOUD 9 FABRICS



**BLACK CHERRY**

**PACIFIC**


**SPRUCE**



**STONE**

**WALNUT**


#### LEXIS FABRICS



**BLACK**

**NEW CHESTNUT**

#### ULTRAFABRICS® ULTRALEATHER™



**PECAN**

**FUDGE**

**GARNET**

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.