# RAVELER

The AmeriGlide Traveler brings exceptional value and easy transporta tion together in a compact, lightweight package. The stylish Traveler has a weight capacity of 275 lbs. and disassembles quickly and easily.



AmeriGlide

## RAVELER

#### Specifications

•	
Max. weight capacity	275 lbs.
Maximum speed <sup>1</sup>	Up to 4 mph
Ground clearance <sup>2</sup>	2.5"
Turning radius <sup>2</sup>	32.5"
Overall length <sup>2</sup>	37"
Overall width <sup>2,6</sup>	19.5"
Tires	
Front	2" x 8" solid
Rear	2.5" x 8" solid
Anti-tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded plastic
Drivetrain	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 2 amp
Per-charge range 1,5	Up to 8 miles (w/ 2 x 12 AH)
Battery requirements <sup>3</sup>	(2) 12V deep-cycle (standard)
Compatible batteries	12 AH
Weight	
Batteries (2 required)	12 AH (19.5 lbs. battery pack) <sup>4</sup>
Weight w/o battery pack <sup>2</sup>	3-wheel: 72 lbs.
Heaviest piece <sup>2</sup>	3-wheel: 28 lbs. (front section)
Standard Colors	

### <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of plus or minus 10%.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

#### **Features**

Easily disassembles into 5 super lightweight pieces for convenient transport and storage

Sleek sporty style

Modular design for easy serviceability

Front frame-mounted seat post offers maximum stability

PDAC coded K0800 power operated vehicle





<sup>&</sup>lt;sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of plus or minus 3%.

<sup>&</sup>lt;sup>3</sup> AGM or Gel-Cell type required.

<sup>&</sup>lt;sup>4</sup> Battery weight may vary based on manufacturer.

<sup>&</sup>lt;sup>5</sup> This range was obtained with the 4-wheel model.

<sup>&</sup>lt;sup>6</sup> Overall scooter width at seat is 22.3"