



Roscoe Mobility Transport Chair with 12" Rear Wheels

The Roscoe Transport Chair features a black aluminum frame that is lightweight and strong. This transport chair comes complete with 12 inch rear wheels and companion loop style brakes that are convenient and easy to use. Black padded nylon upholstery is lightweight and easy to clean. The fixed, full-length, padded arms provide additional comfort.

Product Features

- Black aluminum frame is lightweight and strong
- Companion activated loop style brakes
- Black padded nylon upholstery is lightweight and easy to clean
- Permanent rear axle
- Back release hinge allows back to fold down for easy storage and transport
- Fixed, full-length, padded arms provide additional comfort
- Comes with seat belt for added safety
- Swing-away, composite footrests with heel loops come standard
- 12" flat free pneumatic rear wheels
- 8" x 1" solid rubber front casters
- Easy push-to-lock rear wheel locks
- Limited Lifetime Warranty

Product Specifications

- Product Weight: 22 lbs.
- Overall Width with Arm: 23½"
- Arm Style: Fixed Full Length
- Overall Depth: 24½" (without front rigging)
- Seat Width: 19"
- Seat Depth: 16"
- Seat-to-Floor Height¹: 19¼"
- Back Height: 16"
- Back Style: Folding
- Weight Capacity: 250 lbs.
- HCPCS Code: E1038*

Item #	Finish	Seat Width	Armrests	Front Riggings
KT1912B	Black	19"	Fixed full-length arms	Swing-away footrests

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.
¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com