

Jazzy[®]

SPORT 2

The Jazzy[®] Sport 2 delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Jazzy Sport 2 stability and maneuverability for solid performance indoors and out.

Jazzy Sport[®] 2
shown with Jazzy[®] Armor[™] (standard)



AmeriGlide
Accessibility
Solutions

Specifications

Weight capacity	300 lbs.
Maximum speed	Up to 3.7 mph
Ground clearance	1.5"
Turning radius	24.75"
Overall length	35" without front riggings
Overall width	23"
Drive wheels	9"
Front wheels	3" solid casters
Rear wheels	6" solid casters
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 8.95 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	18 lbs. each
Base weight	68.5 lbs.
Standard seat weight	43.5 lbs. (medium-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Features

- In-line, front-wheel drive technology
- 40 amp, PG GC3 controller
- Easy access to batteries from rear
- Simple main frame design for easy serviceability
- Comfort high-back seat with headrest
- Jazzy® Armor includes shroud and controller guards to protect from daily wear and tear

Options

- High-back seat with solid seat pan
- Cane/crutch holder
- Walker holder
- Cup holder
- Weather cover
- Deluxe Jazzy® Armor™ Package (call for availability)
- Rear basket
- Oxygen tank holder

Recommended Lifts

- Backpacker® Plus
- Outlander
- Backpacker® MV



Color Choices



Jazzy Red