

Uplift Commode Assist:

The lifting commode chair that gives independence and privacy when it's needed most.

Uplift Commode Assist is a self-powered, lifting commode chair that helps users gently seat and raise themselves. It provides safe, controlled support and lifts up to 80% of a user's weight.

Versatile and sturdy, the Commode Assist can be used standalone or over the toilet. For safety, it automatically locks into position when fully depressed.

Benefits

- *Helps keep people active and independent.*
- *Relieves stress and strain on joints.*
- *Lifts only as needed, allowing users to maintain muscular strength.*

Who can benefit

- *People who have lost sufficient upper or lower body strength to get up unassisted.*
- *Those who have lost mobility as a result of:*
 - *Muscular dystrophy*
 - *Lou Gehrig's disease (ALS)*
 - *Post-operative conditions*
 - *Degenerative joint disease*
 - *Arthritis*
 - *Parkinson's disease*



Bringing quality to life.



The Commode Assist eases you gently into a seated position and automatically locks in place. To stand, simply press down on the seat release lever (see below).



Using both hands to support yourself, simply shift forward and begin to stand.



The lifting mechanism activates automatically and provides a gentle and stable lift up from the commode.



Easy to clean snap-on seat.

Removable back rest allows for easy positioning over toilet.

Self-powered lifting mechanism releases slowly as user stands.

Seat release lever & auto-lock hook - commode assist locks when depressed.

Cushioning pads enhance comfort.

Non-skid rubber tips for added safety.

Extendable legs with five height settings.

Bucket slides out for easy cleaning. Includes splash guard and cover.

Steel frame is white powder coated to protect against moisture.

Commode Assist Specifications

| | | |
|-----------------|----------|----------------|
| Base Dimensions | Depth: | 21"/53 cm |
| | Width: | 26"/66 cm |
| Seat Dimensions | Depth: | 16.75"/43 cm |
| | Width: | 16.5"/42 cm |
| Seat Height | Minimum: | 17.75"/45 cm |
| | Maximum: | 23.5"/60 cm |
| Lift | | 9"/23 cm @ 20° |
| Weight | | 21 lb/9.5 kg |



To release the auto-lock, simply push down on the seat release lever. Now both hands are free to use for support when rising (see above sequence).



| Model No. | Description | Weight Range |
|-----------|-----------------------|--------------------------|
| CA 200* | Lifting commode chair | 80-300+ lb 35-136+ kg |

*some assembly required



Bringing quality to life.

Call Toll Free: 1-800-387-0896
 or (902) 422-0804
 Fax Toll Free: 1-888-355-6660
 or (902) 422-0798
 Email: info@up-lift.com

Uplift Technologies Inc.
 11 Morris Drive, Suite 125
 Dartmouth, NS, Canada
 B3B 1M2
 www.up-lift.com