



KM-5000 Series

*Reclining, high strength, lightweight
Aluminum, detachable-arm
wheelchairs*

- Folding 6061 T-6 Aircraft-grade aluminum frame in black powder coat finish
- Medicare Code: K0004+E1226+E0966
- Backrest Reclines from 90°~163° with anti-tippers
- Padded detachable armrests with side panels
- Swing-away detachable elevating leg rests
- Adjustable length leg supports & footplates
- Adjustable-height pillow
- 24" x 1" mag wheels with flat free tires
- Flat free front casters
- Adjustable toggle wheel locks.
- High-grade flame retardant, breathable nylon upholstery
- Limited lifetime warranty on frames and

cross braces, one year on parts*.

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-5000-16	16	18	19.5	33 (18.5+14.5)	50.5	25	36	250
KM-5000-18	18	18	19.5	33 (18.5+14.5)	50.5	27	36	250
KM-5000-20	20	18	19.5	33 (18.5+14.5)	50.5	29	38	300
KM-5000-22	22	18	19.5	33 (18.5+14.5)	50.5	31	38	350

Model No.	Description
KM-5000-16	16", aluminum, 180° reclining back, detachable arms, elevating leg rests, anti-tippers, adj. height pillow
KM-5000-18	18", aluminum, 180° reclining back, detachable arms, elevating leg rests, anti-tippers, adj. height pillow
KM-5000-20	20", aluminum, 180° reclining back, detachable arms, elevating leg rests, anti-tippers, adj. height pillow
KM-5000-22	22", aluminum, 180° reclining back, detachable arms, elevating leg rests, anti-tippers, adj. height pillow

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Leg Rests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Quick Release Axle • Spoke Wheels • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Tel.626-581-2235 Fax.626-581-2335