

ORFIT<sup>®</sup>

EASE



# ORFIT® EASE

## Orfit Industries surpasses the classic plastic and rubber (-like) splinting materials

This rigid splinting material is designed for the fabrication of larger splints. It offers users high resistance to stretch, combined with the highest rigidity and durability available in a splinting material. The performance characteristics of ORFIT® EASE make it the material of choice for an unsurpassed range of splints. Its versatility is unmatched.

Extremely easy to form, ORFIT® EASE becomes extremely hard after cooling down. It has no elastic behaviour and it will not recoil during moulding.

The product's ease of moulding is the prime reason for its widespread use by so many professionals.

Our analysis of the market confirms that traditional splinting materials are still in demand and widely used. As a result of this assessment, Orfit Industries has developed the "New Classics" to include in our growing family of high quality splinting materials.

The "New Classics" from Orfit Industries include ORFIT® ORTHO, ORFIT® DRAPE NS and ORFIT® EASE...

They provide you with the "Best of all Worlds". You get an improved and superior surface finish, while maintaining the "feel" that you are accustomed to. The "New Classics" do not contain rubber and are not sensitive to UV light.

Orfit Industries is now able to offer you a complete selection of "memory" and "non memory" splinting materials, with the added benefit of Orfit® quality that is recognized and respected worldwide.

### The "New Classics", Made Better by Orfit Industries

ORFIT® EASE			
Thickness	Perforation	Size	
		450 x 600 mm	600 x 900 mm
1,6 mm	Non perforated	1232.1	
	Macro	1232.6	
3,2 mm	Non perforated	1234.1	
	Macro	1234.6	
	Non perforated		1254.1
	Macro		1254.6



ORFIT Industries  
Vosveld 9A  
B-2110 Wijnegem  
Belgium

Phone: (+32) (0)3.326.20.26  
Fax: (+32) (0)3.326.14.15  
www.orfit.com