



KN-920W

Standard steel bariatric

- *High quality steel construction*
- *Detachable desk & full armrests*
- *Swing-away footrests*
- *24" x 1" mag wheels w/ solid tires*
- *Double padded vinyl upholstery*
- *Frame: Chrome*
- *Weight Capacity: 300 lbs.*

*

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/o riggings (lb.)	weight capacity (lb.)
KN-920W	20"	16"/18"	21"	18"	36"	28"	N/A	300
KN-920WE	20"	16"/18"	21"	18"	36"	28"	N/A	300

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Leg Straps
- Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335