



LT-800B Series

*Standard lightweight detachable arm
Wheelchairs*

- Medicare code k0003
- 18" Standard Seat width
- Fixed padded full armrests
- Swing away footrests
- 24" Flat free polyurethane rear tires
- 8" Flat free front casters
- Double padded black nylon upholstery
- Adjustable seat height
- W/Heel loop
- Colors: Black only
- Wheelchair Weight: 35 lbs (w/o footrests)
- Weight Capacity: 250 lbs

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/o riggings (lb.)	weight capacity (lb.)
LT-800B	18	16	20	17	36	26.5	35	250
LT-800NB	16	16	20	17	36	24.5	36	250

Model No.	Description
LT-800B	18", Padded fixed armrests, swing-away footrests
LT-800NB	16", Same as above

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Leg Straps
- Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335