

# S-Ergo Lite











## Ultra Lightweight Ergonomic Transporter



- LIGHTEST TRANSPORTER IN THE MARKET!
- S-Style Ergonomic Seat
- High strength, ultra-lightweight only 17lbs. (w/ footrests)
- Built in Silver Anti-bacterial 1/4" cushion
- Folding backrest/folding seat for transporting in vehicle or travel
- Fixed armrests w/ wider concaved armpads
- Fixed footrests w/ extra wide footplate
- Push bar brake lever for attendant, also lockable
- Pouch for carrying small items attached to upholstery behind backrest & small carry pouch on each armrest
- 14" Rear polyurethane, high tread, flat free wheels
- 7" x 1" Polyurethane front caster
- 16" x 17" OR 18" x 17" Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black/Silver
- Weight Capacity of 250 lbs

### Features for S-Ergo Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slipage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/concaved armpad, providing comfort & support

items	F/R wheels											Max. User Weight
Width												
S-Ergo Lite ~ 16"x17"	7"/14"	16"	17"	23"	12"	8"	35 3/4"	18 1/2"	16"	36"	17 lbs.	250 lbs.
S-Ergo Lite ~ 18"x17"	7"/14"	18"	17"	25"	12"	8"	35 3/4"	18 1/2"	16"	36"	17 lbs.	250 lbs.

*\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam, UV Resistant Canopy, Rigidifying Push-Bar, Foldable/Adjustable Headrest, Frog Leg Suspension*



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue, City of Industry, CA 91748-1418

Order. 1-800-805-2762 / Tel.626-581-2235 / Fax.626-581-2335