

Invacare® P9000™ XDT Power Wheelchair



Yes, you can.®

P9000™ XDT





Yes, you can.®

The Invacare® P9000™ XDT is a traditional rear-wheel drive chair that maneuvers easily and negotiates obstacles both inside and outside the home. It offers all the features consumers want at an affordable price. And it easily disassembles and folds for transport.

Features

- Rear-wheel drive and compact size provide optimum maneuverability and accessibility
- Invacare® MK₅™ NX electronics featuring the CSPJ joystick
- Folds to narrow 13" for transport
- Adjustable-height back provides comfort and support
- Powerful 350-watt, high-torque motors ensure dependable performance
- Telescoping front riggings allow optimal foot placement
- Wheel locks for added safety
- 250 lb. user weight capacity
- Available in two seat sizes: 18" W x 16" D, 18" W x 18" D



Specifications (All weight/measurements are approximate.)

Dimensions	
Length	47"
Width (excluding joystick)	24.5"
Height	38"
Seat-to-floor height	19.5"
Seat Dimensions	
Width	18"
Depth	16", 18"
Back height	17", 18", 19"
Seat angle	3°
Weight capacity	250 lb.
Wheel & Caster Sizes	
Front/Rear casters	8" x 1.75"
Drive wheels	12.5" x 2.25"
Maximum Speed	
	4.0 mph
Incline Capability	
	9°
Turning Radius	
Front (footrest in place)	36"
Rear (anti-tippers in place)	16"
Weight	
*with UI batteries	142 lb.
without batteries and front riggings	82 lb.
Ground Clearance	
	3.25"
Limited Warranties	
Base frame	36 months
Seat frame	36 months
Electronics	12 months
Motor/Gearbox	12 months

*Measured with front riggings and anti-tippers with 16" seat depth.



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2007 Invacare Corporation. All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.

Invacare Corporation
Form No. 03-067
Rev. 10/07