

Invacare® Tracer® EX2 wheelchair sets the benchmark in manual wheelchairs. A true dual axle position allows for repositioning the 24" rear wheels and 8" casters for adult and hemi seat-to-floor heights. Offset permanent arms offer full seat width up to 20" with a permanent arm, hemi-spaced footrest and full range of accessories. Incorporating the design and technology of the Invacare 9000 series, the Tracer EX2 wheelchair has interchangeable components with the 9000 series chairs.



Invacare® Tracer® EX2 Wheelchair

Features

- Durable, low-maintenance, triple chrome-plated, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- Offset permanent arms available in 16", 18" and 20" widths
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Flame-retardant upholstery resists mildew and bacteria, as well as cold-crack and UV damage
- Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

Specifications

Seat sizes	
Widths	16", 18", 20"
Depths	16" standard
Seat-to-floor heights	
	17.5", 19.5"
Back height	
	16" fixed
Casters and wheels	
Casters	8"
Wheels	24"
Front riggings	
Hemi swingaway footrest	17"-21" knee-to-heel
Hemi swingaway padded elevating legrest	17.5"-21" knee-to-heel
Articulating elevating legrest	17"-25" knee-to-heel
(Footplates available in aluminum or composite)	

Arm styles	
	Permanent full-length, padded
	Full-length fixed or adjustable height, padded, removable
	Desk-length fixed or adjustable height, padded, removable
	(All arms are conventional)
Weight capacity	
	250 lb.
Product weight (TREX2/WD86/28)	
	36 lb.
Limited warranty	
	5 years on frame and crossbraces



Offset permanent arms allow for full seat width between arms.