



KN-700T

Standard steel detachable armrest wheelchair

- Medicare code K0001
- Detachable/Padded desk length armrests with alloy side guard
- Swing-away/Height adjustable footrests with heel loop
- Dual-Axle/Seat height adjustable frame
- Push to lock wheel brakes
- 24" Mag wheels with flat free tires
- Extra-durable forks
- 8" Flat free front casters
- High-grade/Padded flame retardant vinyl upholstery
- Carry pocket on backrest
- Seat width: 16", 18"
- Frame Color: Silver vein
- Weight Capacity: 250 lbs

SPECIFICATIONS (subject to change without notice) :

Model No.	Seat Width (in.)	Seat Depth (in.)	Armpad Height (in.)	Seat Height (in.)	Back Height (in.)	Overall Height (in.)	Overall Length (in.)	Folded Width (in.)	Overall Width (in.)	Weight w/o Riggings (lbs.)	Weight Capacity (lbs.)
KN-700NT	16	16	9	17 or 19	16	36	42(31)	11	24	39	250
KN-700T	18	16	9	17 or 19	16	36	42(31)	11	26	39	250

Optional Accessories:

- Elevating Legrests - Anti-Tippers - Brake Extension - IV Pod Attachment - Oxygen Tank Holder - Seat Belt - Leg Straps
- Carry Pouch - Headrest - Backrest Extension



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335