



Yes, you can.®

The Invacare® Veranda™ Wheelchair provides comfort and mobility in an economical package. It comes complete with front riggings packaged in the carton with the chair. The Invacare Veranda combines style and durability at an affordable price.



Features

- Padded nylon upholstery
- Comfortable padded arms
- Stylish powder-coated frame includes front riggings
 - Footrests match chair
 - Legrests are chrome
- Carry pocket on back
- Large step-tube end cap

Specifications

Seat sizes	
Width	18"
Depth	16"
Seat-to-floor height	
Removable arm	19.5"
Permanent Arm	19.5"
Back height	18"
Front riggings	
Hemi swingaway footrest	16.5" - 19.5" knee-to-heel
Hemi swingaway padded elevating legrest	17.5" - 21" knee-to-heel
Arm styles	
Permanent full-length, fixed height arm	
Removable desk-length, fixed height arm	
Weight capacity	250 lb.
Product weight	37 lb.
Limited warranty	1 year



Pocket on back.



Large step-tube end cap.



Step-tube.

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.

© 2007
Invacare Corporation

Form No. 07-064

Invacare® Veranda™ Wheelchair



Veranda™ Wheelchair

Yes, you can.®

