

Burgundy Finish



Blue Finish



Roscoe Mobility Aluminum Transport Chair

The Roscoe Aluminum Transport Chair is lightweight, strong and available in your choice of blue or burgundy finish. The back release hinge allows the back to fold down for easy storage and transport. Weighing less than 20 pounds, it is equipped with comfortable, padded nylon upholstery. Comes complete with seat belt for added safety. Industry compatible swing-away, composite footrests with heel loops come standard.

Product Features

- Aluminum frame is lightweight and strong
- Permanent rear axle
- Black padded nylon upholstery is lightweight and easy to clean
- Back release hinge allows back to fold down for easy storage and transport
- Fixed, full-length, padded arms provide additional comfort
- Comes with seat belt for added safety
- Swing-away, composite footrests with heel loops come standard
- Composite 8" x 1" wheels with solid rubber tires are lightweight and maintenance free
- Easy push-to-lock rear wheel locks
- Limited Lifetime Warranty

Product Specifications

- Product Weight: 19½ lbs.
- Overall Width with Arm: 23½"
- Arm Style: Fixed
- Overall Depth: 27" (without front rigging)
- Seat Width: 19"
- Seat Depth: 16"
- Seat-to-Floor Height¹: 19"
- Back Height: 16"
- Back Style: Folding
- Weight Capacity: 250 lbs.
- HCPCS Code: E1038*

Item #	Finish	Seat Width	Armrests	Front Riggings
KTA1916SA-BG	Burgundy	19"	Fixed full-length arms	Swing-away footrests
KTA1916SA-BL	Blue	19"	Fixed-full length arms	Swing-away footrests

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice. ¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.



To Order:
(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com