MODEL

Jazzy 614HD 2HD-SS ¹	Group 2: Heavy Duty (K0824)
Jazzy 614HD 2HD-C ¹	Group 2: Heavy Duty (K0825)
Weight Capacity	450 lbs.
Turning Radius	22.75"
Overall Width ³	25.25"
Overall Length ³	35.5″
Maximum Speed ²	Up to 4.25 mph
Ground Clearance ³	2.25″
Front Wheels	6" solid casters
Drive Wheels	14" foam filled
Rear Wheels	6" solid casters
Drivetrain	Two-motor, in-line, Mid-Wheel 6® drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Active-Trac [®] (ATX)
High-back Seat Sizes	18"x 18"-20", 20" x 18" - 20", 22" x 20" - 22"
Seat-to-Floor Range	20.5" - 23.5"
Per Charge Range ^{3,5}	Up to 11 miles
Standard Electronics	50A, Dynamic Shark
Component Weights	Base: 122 lbs. (heaviest piece) Seat: 18"- 51.8 lbs.; 20"- 54 lbs. 22"- 57.6 lbs. (comfort high-back) Batteries: 37 lbs. each
Battery Requirements ^{4,6}	(2) 12 volt, NF-22
Battery Charger	5 amp, off-board
Warranty	Frame: Lifetime limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

1 Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

6 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



(US) 800-800-8586 • (Canada) 888-570-1113 pridemobility.com | 🖸 🔂 🐻

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.



POWER AND PRECISION WITH A **TIGHT** TURNING RADIUS



Medical Device¹

AMERICA'S FAVORITE POWER CHAIR gojazzy.com | 🖸 f 😏 💿

HEAVY-DUTY PERFORMANCE



FDA Class II Medical Device¹

Mid-Wheel 6[®] Technology allows for maximum stability



Active-Trac[®] ATX Suspension



Easy-access freewheel levers



50A, dynamic shark controller

Power and precision combine into one exceptional power chair! The Jazzy[®] 614 HD¹ features outstanding suspension and a tight turning radius. Large drive wheels and six-inch casters tackle uneven terrain. no matter where your day takes you.

- Universal, fixed, in-line joystick mount
- Lap belt

Available Colors



Candy Viper Blue Apple Red

Seating Sizes

High-back, black vinyl seat up to 24"W x 24"D

Accessories

- Cane/crutch holder
- Cup holder

Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away legrests •
- Swing-away joystick

- Flag holder
- Lap belt (50", 60", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover